


Product Datasheet

Human ZNF804A (Zinc finger protein 804A) ELISA Kit (orb1100969)

Description	Human ZNF804A (Zinc finger protein 804A) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	ZNF804A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7Z570
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
ZNF804A
(Zinc
finger
protein
804A)...