


Product Datasheet

Human ROS1 (Proto-oncogene tyrosine-protein kinase ROS)ELISA Kit (orb1100970)

Description	Human ROS1 (Proto-oncogene tyrosine-protein kinase ROS)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	ROS1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P08922
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
ROS1
(Proto-
oncogene
tyrosine-
prot...