


Product Datasheet

Human SIRT3(NAD-dependent ADP-riboseyltransferase sirtuin-4) ELISA Kit (orb1100975)

Description	Human SIRT3(NAD-dependent ADP-riboseyltransferase sirtuin-4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SIRT3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O75094
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
SIRT3(NAD-
dependent
ADP-
riboseyltra...