


Product Datasheet

Human SLIT3(Slit homolog 3 protein) ELISA Kit (orb1100976)

Description	Human SLIT3(Slit homolog 3 protein) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	SLIT3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O75094
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
SLIT3(Slit
homolog
3 protein)
ELIS...