


Product Datasheet

Human CACNA2D1 (calcium channel, voltage-dependent, alpha 2/delta subunit 1) ELISA Kit (orb1100982)

Description	Human CACNA2D1 (calcium channel, voltage-dependent, alpha 2/delta subunit 1) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CACNA2D1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P54289
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human CACNA2D1 (calcium channel, voltage-dependent, alpha 2/delta subunit 1)