


Product Datasheet

Human ZNRD1 (DNA-directed RNA polymerase I subunit RPA12) ELISA Kit (orb1100989)

Description	Human ZNRD1 (DNA-directed RNA polymerase I subunit RPA12) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	ZNRD1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9P1U0
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
ZNRD1 (DNA-
directed RNA
polymerase...