


Product Datasheet

Human TYRP1 (5,6-dihydroxyindole-2-carboxylic acid oxidase) ELISA Kit (orb1100991)

Description	Human TYRP1 (5,6-dihydroxyindole-2-carboxylic acid oxidase) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	TYRP1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P17643
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human TYRP1
(5, 6-
dihydroxyindole-
2-carb...