


Product Datasheet

Human TBK1 (Serine/threonine-protein kinase TBK1) ELISA Kit (orb1100992)

Description	Human TBK1 (Serine/threonine-protein kinase TBK1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	TBK1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UHD2
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human TBK1
(Serine/threonine-protein kin...