


Product Datasheet

Human FUNDC1 (FUN14 domain-containing protein 1) ELISA Kit (orb1101010)

Description	Human FUNDC1 (FUN14 domain-containing protein 1) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	FUNDC1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8IVP5
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
FUNDC1
(FUN14
domain-
containing
pr...