


Product Datasheet

Human C3AR1 (C3a anaphylatoxin chemotactic receptor) ELISA Kit (orb1101013)

Description	Human C3AR1 (C3a anaphylatoxin chemotactic receptor) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	C3AR1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q16581
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
C3AR1 (C3a
anaphylatoxin
chemotact...