


Product Datasheet

Human CBR1 (Carbonyl reductase [NADPH] 1) ELISA KIT (orb1101015)

Description	Human CBR1 (Carbonyl reductase [NADPH] 1) ELISA KIT
Reactivity	Human
Range	0.156-10ng/ml
Target	CBR1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P16152
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
CBR1
(Carbonyl
reductase
[NADPH]
1...