


Product Datasheet

Human AKT1S1 (Proline-rich AKT1 substrate 1) ELISA Kit (orb1101016)

Description	Human AKT1S1 (Proline-rich AKT1 substrate 1) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	AKT1S1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96B36
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
AKT1S1
(Proline-
rich AKT1
substrat...