


Product Datasheet

Human β -ENaC(Amiloride-sensitive sodium channel subunit beta) ELISA Kit (orb1101019)

Description	Human β -ENaC(Amiloride-sensitive sodium channel subunit beta) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	β -ENaC
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P51168
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human β -ENaC(Amiloride-sensitive sodium...)