


Product Datasheet

Human ZNF746(Zinc finger Protein 746) ELISA Kit (orb1101020)

Description	Human ZNF746(Zinc finger Protein 746) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	ZNF746
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q6NUN9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
ZNF746(Zinc
finger
Protein 746)
EL...