


Product Datasheet

Human VAPB(Vesicle-associated membrane protein-associated protein B/C) ELISA Kit (orb1101023)

Description	Human VAPB(Vesicle-associated membrane protein-associated protein B/C) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	VAPB
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O95292
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
VAPB(Vesicle-
associated
membrane
p...