


Product Datasheet

Human TERF2 (Telomeric repeat-binding factor 2) ELISA Kit (orb1101030)

Description	Human TERF2 (Telomeric repeat-binding factor 2) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	TERF2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15554
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
TERF2(Telomeric
repeat-binding
fac...