


Product Datasheet

Human SNAP23(Synaptosomal-associated protein 23) ELISA Kit (orb1101037)

Description	Human SNAP23(Synaptosomal-associated protein 23) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	SNAP23
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O00161
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
SNAP23(Synaptosomal-associated pro...