


Product Datasheet

Human SIRT2(NAD-dependent protein deacetylase sirtuin-2) ELISA Kit (orb1101039)

Description	Human SIRT2(NAD-dependent protein deacetylase sirtuin-2) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	SIRT2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8IXJ6
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
SIRT2(NAD-
dependent
protein
deacet...