


Product Datasheet

Human SARS2(Serine-tRNA ligase, mitochondrial) ELISA Kit (orb1101040)

Description	Human SARS2(Serine-tRNA ligase, mitochondrial) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	SARS2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9NP81
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
SARS2(Serine-
tRNA ligase,
mitochon...