


Product Datasheet

Human PSMA6 (Proteasome subunit alpha type-6)ELISA Kit (orb1101050)

Description	Human PSMA6 (Proteasome subunit alpha type-6)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	PSMA6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P60900
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
PSMA6
(Proteasome
subunit
alpha ty...