


Product Datasheet

Human P2RY6 (P2Y purinoceptor 6) ELISA Kit (orb1101056)

Description	Human P2RY6 (P2Y purinoceptor 6) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	P2RY6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15077
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
P2RY6 (P2Y
purinoceptor
6) ELISA K...