


Product Datasheet

Human OSGEPL1(Probable tRNA N6-adenosine threonylcarbamoyltransferase, mitochondrial) ELISA Kit (orb1101057)

Description	Human OSGEPL1(Probable tRNA N6-adenosine threonylcarbamoyltransferase, mitochondrial) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	OSGEPL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H4B0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
OSGEPL1(Probable
tRNA N6-
adenosine...