


Product Datasheet

Human OSBP(Oxysterol-binding protein 1) ELISA Kit (orb1101059)

Description	Human OSBP(Oxysterol-binding protein 1) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	OSBP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P22059
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
OSBP(Oxysterol-
binding protein
1) ...