


Product Datasheet

Human OGDHL (2-Oxoglutarate Dehydrogenase-Like) ELISA Kit (orb1101061)

Description	Human OGDHL (2-Oxoglutarate Dehydrogenase-Like) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	OGDHL
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9ULD0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human OGDHL
(2-Oxoglutarate
Dehydrogenas...