


Product Datasheet

Human NSFL1C(NSFL1 cofactor p47) ELISA Kit (orb1101063)

Description	Human NSFL1C(NSFL1 cofactor p47) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	NSFL1C
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UNZ2
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
NSFL1C(NSFL1
cofactor p47)
ELISA K...