


Product Datasheet

Human MGAT2 (N-glycosyl-oligosaccharide-glycoprotein N-acetylglucosaminyltransferase)ELISA Kit (orb1101068)

Description	Human MGAT2 (N-glycosyl-oligosaccharide-glycoprotein N-acetylglucosaminyltransferase)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	MGAT2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q10469
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human MGAT2
(N-glycosyl-oligosaccharide-
...)