


Product Datasheet

Canine MMP-9(Matrix Metalloproteinase 9) ELISA Kit (orb1101078)

Description	Canine MMP-9(Matrix Metalloproteinase 9) ELISA Kit
Reactivity	Canine
Range	0.156-10ng/ml
Target	MMP-9
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O18733
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Canine MMP-9(Matrix Metalloproteinase 9)...