


Product Datasheet

Human HS-IL-6(High sensitive Interleukin 6) Accquant ELISA Kit (orb1101090)

Description	Human HS-IL-6(High sensitive Interleukin 6) Accquant ELISA Kit
Reactivity	Human
Range	0.078-5pg/ml
Target	IL-6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P05231
Sensitivity	0.046pg/ml
Expiration Date	Please enquire.



Human HS-IL-6(High sensitive Interleukin...)