


Product Datasheet

Human HS-IL-15 (High sensitive Interleukin 15) Accquant ELISA Kit (orb1101092)

Description	Human HS-IL-15 (High sensitive Interleukin 15) Accquant ELISA Kit
Reactivity	Human
Range	0.31-20pg/ml
Target	IL-15
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.875pg/ml
Expiration Date	Please enquire.



Human
HS-IL-15
(High
sensitive
Interleuk...