


Product Datasheet

KIR2DL1 antibody (orb1101101)

Description	Mouse monoclonal antibody to KIR2DL1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS
Immunogen	The immunogen consisted of KIR2DL1 ectodomain fused to the hIgG1 Fc fragment and was produced from the supernatant of HEK 293T cells
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



Dot plot graph showing CD56+ hPBMCs stained...