


Product Datasheet

m123/IE1 (MCMV) antibody (orb1101111)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to m123/IE1 (MCMV) |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FACS, IF, IHC, WB |
| Immunogen | The immunogen consisted of MCMV virions |
| Form/Appearance | Liquid |
| Concentration | 1 mg/mL |
| Storage | Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | Affinity purified |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of mouse li...

Image not found or type unknown

Immunohistochemical staining of mouse li...