

Product Datasheet


Human Activin RIIA / ACVR2A Protein (orb1101193)

Description

Activin receptor type-2A (ACVR2A) is also known as

Activin Receptor Type 2A (ACVR2A)
1.1 µg of Activin Receptor Type 2A (ACVR2A) per well

Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Activin RIIA
Preservatives	51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 39.9 kDa. The protein migrates as 50-60 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human ACVR2A, Fc Tag (orb1101193) is expressed from human 293 cells (HEK293). It contains AA Ala 20 - Pro 135 (Accession # P27037-1).
MW	39.9 kDa
NCBI	NP_001607.1
Expiration Date	6 months from date of receipt.



Immobilized
Human
Activin A,
premium
gra...