

Product Datasheet

Marijuana antibody (orb11017)

Description	Marijuana antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Cannabinol
Target	Marijuana(2B1)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	0.31043 kDa
CAS Number	521-35-7
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.