

Product Datasheet

IL8 antibody (HRP) (orb110251)

Description

Rabbit polyclonal antibody to IL8 (HRP)

Species/Host

Rabbit

Reactivity

Human, Rabbit

Conjugation

HRP

Tested

ELISA, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from human IL-8

Target

CXCL8

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Application notes

orb110251 is the HRP conjugated form of orb10915

Isotype

IgG

Clonality

Polyclonal

MW

8/9kDa

Uniprot ID
P10145
Dilution Range

IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.