



Product Datasheet

Mfn1 antibody (orb11040)


Description	Mfn1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Mfn 1 (651-741/741 aa)
Target	Mitofusin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	Q8IWA4
Dilution Range	Flow-Cyt=1µg/Test, WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart (la...



Western blot analysis of mouse heart cel...



Flow cytometric analysis of Rsc96 cell u...