

Product Datasheet

SLC26A4 antibody (HRP) (orb110416)

Description

Rabbit polyclonal to SLC26A4 (HRP).

Species/Host

Rabbit

Reactivity

Human

Conjugation

HRP

Tested

ELISA, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Solute carrier family 26 member 4

Target

SLC26A4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Application notes

orb110416 is the HRP conjugated form of orb101538

Isotype

IgG

Clonality

Polyclonal

MW

93kDa

Uniprot ID
[O43511](https://www.uniprot.org/entry/O43511)
Dilution Range

WB=1:500-2000,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.