



Product Datasheet

MIP1 alpha antibody (orb11052)


Description	MIP1 alpha antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MIP-1 alpha (41-92/92 aa)
Target	CCL3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8 kDa
Uniprot ID	P10147
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of rat exfo...



IHC-P of mouse brain tissue (1:200 dilut...



paraffin-embedded mouse brain tissue usi...