



Product Datasheet

MMP1 antibody (orb11056)


Description	MMP1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse MMP1 (27 kDa interstitial collagenase) (401-464/464 aa)
Target	MMP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27/41/54 kDa
Uniprot ID	P03956
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver lys...



Immunohistochemical analysis of formalin...



IHC-P of mouse endometrium tissue (MMP1 ...)