

Product Datasheet

PPAR-gamma antibody (orb1105899)

Description	Rabbit polyclonal antibody to PPAR-gamma
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PPAR-gamma
Preservatives	PBS
Form/Appearance	Liquid
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.