

Product Datasheet

Cleaved-Caspase 3 antibody (orb1105907)

Description	Rabbit polyclonal antibody to Cleaved-Caspase 3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Target	Cleaved-Caspase 3
Preservatives	PBS
Form/Appearance	Liquid
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.