

Product Datasheet

SARS-CoV-2 Nucleoprotein antibody [Out of stock] (orb1106161)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Nucleoprotein
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	LP-13473, synthetic peptide
Preservatives	PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA, WB
Expiration Date	12 months from date of receipt.