

## Product Datasheet

### MMP3 antibody (orb11062)

## Description

MMP3 antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Porcine, Rat

### Conjugation

Unconjugated

### Tested Applications

ICC, IF, IHC-P, WB

### Immunogen

KLH conjugated synthetic peptide derived from human MMP3. Please contact us for the exact immunogen sequence. The peptide is available as orb12895.

### Target

MMP3

### Form/Appearance

10 mM PBS, 0.02% sodium azide

### Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

### Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Polyclonal antibodies are purified by peptide affinity chromatography

### MW

54 kDa

### Uniprot ID

[P08254](#), [P28862](#), [P03957](#)

### NCBI

[NM\\_002422.3](#), [NP\\_002413.1](#), [NM\\_002422.31](#), [NP\\_002413.11](#)

### Entrez


[4314](#)

### Dilution Range


IF/ICC: 1:50-400, IHC-P: 1:50-400, WB: 1:200-2000

### Expiration Date


12 months from date of receipt.



WB analysis of rat kidney (lane 1), mous...



IHC-P image of mouse kidney tissue using...



IHC-P staining of rat spinal cord tissue...