

Product Datasheet

Txn Termination Factor antibody [Out of stock] (orb1106236)

Description	Rabbit polyclonal antibody to Txn Termination Factor
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	A synthesized peptide derived from human TTF1.
Target	TTF1
Preservatives	PBS, pH 7.4, 0.02% Sodium Azide, 50% Glycerol
Concentration	1 mg/ml
Storage	Upon receipt, store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Entrez	7270
Dilution Range	IHC-P (10 µg/ml), Peptide-ELISA (1:20000 - 1:40000), WB (1:1000 - 1:3000)
Expiration Date	12 months from date of receipt.