

Product Datasheet

KLH antibody (orb11067)

Description	KLH antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH protein (please contact us for details)
Target	KLH
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Dilution Range	WB:1:100-500, IHC-P:1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.