



Product Datasheet

Mesothelin antibody (orb11082)


Description	Mesothelin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Mesothelin (101-200/630 aa)
Target	Mesothelin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	64 kDa
Uniprot ID	Q13421
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human pancreatic cancer using M...



Immunohistochemical analysis of formalin...