

Product Datasheet

Histamine H4 Receptor antibody [Out of stock] (orb1108301)

Description	Rabbit polyclonal antibody to Histamine H4 Receptor
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic 21 amino acid peptide from 3rd cytoplasmic domain of human HRH4 / Histamine H4 Receptor. Percent identity with other species by BLAST analysis: Human, Gibbon (100%); Marmoset (90%); Monkey (81%).
Target	HRH4
Preservatives	PBS, 0.1% Sodium Azide
Concentration	0.7 mg/ml
Storage	Aliquot and store undiluted at -20°C or below for up to 1 year. Can be stored undiluted at 4°C for up to 1 month. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Entrez	59340
Dilution Range	IHC, IHC-P (14 µg/ml)
Expiration Date	12 months from date of receipt.