



Product Datasheet

Melatonin Receptor 1A antibody (orb11085)


Description	Melatonin Receptor 1A antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Melatonin Receptor 1A. Please contact us for the exact immunogen sequence. The peptide is available as orb374736.
Target	Melatonin Receptor 1A
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	38 kDa
Uniprot ID	B7X945, P48039, Q61184
NCBI	31, 005958, 005949, 11
Entrez	4543
Dilution Range	WB: 2 µg/ml, IHC-P:1:500, IF/ICC: 1:300
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of mouse hippocampus tissue using ...



IHC-P of mouse hippocampus tissue (orb11...