

Product Datasheet

SARS-CoV-2 Spike S1 Indian B.1.617.2 (delta) (orb1108512)

Description	SARS-CoV-2 Spike S1 Indian B.1.617.2 (delta)
Preservatives	PBS pH 7.4
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze- thaw cycles. Shipping on dry ice.
Note	For research use only
Source	CHO-based cell line (expressed by QMCF Technology)
MW	75.70 kDa, Protein contains amino acids 14-681, mutations T19R, G142D, E156G, del_F157-R158 L452R, T478K, D614G, P681R, two extra amino acids (AS) in N-terminus and His-6 tag at C-terminus and GSG linker between protein and tag.
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.