

Product Datasheet

SARS-CoV-2 Spike S1 Indian B.1.617.1 (kappa) (orb1108513)

Description	SARS-CoV-2 Spike S1 Indian B.1.617.1 (kappa)
Preservatives	PBS pH 7.4
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles. Shipping on dry ice.
Note	For research use only
Source	CHO-based cell line (expressed by QMCF Technology)
MW	79.99 kDa, Protein contains amino acids 14-681, mutations G142D, E154K, L452R, E484Q, D614G, P681R, two extra amino acids (AS) in N-terminus and His-6 tag at C- terminus and GSG linker between protein and tag.
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.