

Product Datasheet

SARS-CoV-2 Spike S1 South Africa VOC 501.V2 (beta) (orb1108516)

Description	SARS-CoV-2 Spike S1 South Africa VOC 501.V2 (beta)
Preservatives	PBS pH 7.4
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles. Shipping on dry ice.
Note	For research use only
Source	CHO-based cell line (expressed by QMCF Technology)
MW	75.454 kDa, protein contains amino acids 14-681, mutations L18F, D80A, D215G, LAL 242-244 del, R246I, K417N, E484K, N501Y, D614G, two extra amino acids (AS) in N-terminus and His-6 tag at C-terminus and GSG linker between protein and tag.
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.